



Gene Check, Inc.
www.GeneCheck.com

**CERVID
SUBMISSION FORM**

1175 58th Ave., Suite 100
Greeley, CO 80634
970-472-9951
genecheck@genecheck.com

OWNER NAME (REQUIRED):	VETERINARIAN (IF APPLICABLE):
FARM NAME (OPTIONAL):	Address:
MAILING ADDRESS:	City, State Zip:
CITY, STATE, ZIP:	VET PHONE: Vet Email:
PHONE:	EMAIL:

<i>DEER DNA Test</i>	<i>Cost</i>
Codon 96 (GG, GS, SS results)	\$18
<i>ELK DNA Test</i>	<i>Cost</i>
Codon 132 (MM, ML, LL results)	\$18

Note: there is a \$2 additional charge per sample per test for hair samples.

PAYMENT METHOD: Check Enclosed Charge my Credit Card (Visa, MasterCard, Discover)

Card #: _____

Exp date: _____ 3 digit security code on back of card _____

Signature: _____

Printed Name: _____

DNA SAMPLE METHODS & SHIPPING INFORMATION

BLOOD TUBES- Collect 2-3 cc's of whole blood in a **"purple top" EDTA tube**
(Be sure to rock the tube back and forth immediately after collection).
Shipping tubes: Ship tubes in well padded box.
In hot months ship via 1-3 day delivery with ice. In cooler months, ship via 1-5 day delivery without ice.

DNA EAR PUNCH systems such as Datamars, Hauptner-Typifix & Allflex (Call the office to order supplies).
These methods all have a preservation agent in the collection tube. After collection the tissue can be stored for several years at room temperature. Do not refrigerate or ice these types of tissue sample.
SHIPPING DNA TISSUE PUNCHES: Ship in a padded envelope or box via regular/ ground delivery.

FROZEN TISSUE:
The tissue can be taken from almost anywhere on the animal, but is usually an ear or tail sample. The sample must be a minimum of 1/4" and **must be frozen very soon after collection.**
If the animal is deceased, any non decayed sample is acceptable.
SHIPPING FROZEN TISSUE: Ship over night on ice.

BLOOD CARDS
We are set up to use FTA Whatman blood cards. If you choose to submit another brand of blood card, we will do our very best to get a result for you, but we can not guarantee that we will be successful. We use the Whatman protocol and reagents and this method is not ideal for other blood card brands. If we are unable to get a result, you would need to resubmit the sample and pay to have the sample tested again. If you use Whatman cards and we do not get a result, you can resubmit at no additional charge.
SHIPPING BLOOD CARDS: Blood cards can be shipped via regular/ ground delivery.
Rush testing is not available for blood cards.

HAIR SAMPLES
At least 30-40 hairs with follicles taped down to paper in a bunch. Do not tape over the follicle. Please make sure hair is clean, free of debris and is grouped together (not loose).
There is a \$2 additional charge per sample for hair samples.

*** Place an "X" under the desired test column for each sample ***

Sample Information: (Please Print Clearly)		Deer Codon 96	Elk Codon 132	*RUSH \$10 per test See note at bottom of page	Notes or Other ID#s (will not be printed on certificate)
Deer/Elk ID# (This # will be printed on your certificate)					
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					

* **RUSH TESTING** 24-48 hour RUSH option is available for DNA testing only (weekends are not included in the 24-48 hours). Cost is \$10 extra per test. Rush testing is not available for blood cards.

Check the "RUSH" box next to each sample that you would like to have rushed. If there is more than one test requested for that sample, indicate in the notes section if all tests should be rushed for that sample, or only some should be rushed.